

Granulometer 920 number 353

11/30/2000 16:56:49

File name : C:\GRANULOG\OP36006.MES

Version : 3.30

Ref. CILAS : 4561C9300197F003C38D1

Sample : OP 360 DEL 18/11/2000

Liquid : EMULSIONE H2O E DETERG. CARECLEAN XHP/F

Ultrasonic mixer : 20 s. / Dispersing agent : ACQUA

Comment : LAVATRICE OP 360

User name : W.M.T. s.r.l. - FR

Plant :

Place :

Concentration : 392

Median size : 3.32  $\mu$

Diameter at 10.0 % : 0.78  $\mu$

Diameter at 90.0 % : 15.95  $\mu$

Cumul. at 100.00  $\mu$  : 100.00 %

Number of measur. : 20

Number of cleaning : 4

automatic dilution : No


US mixer / measur. : Yes

Option : -----

To top up : No

Results : Weight distribution / Undersize

D	0.30	0.60	1.00	1.40	1.70	2.00	2.60	3.20	4.00	5.00
C%	0.00	3.35	16.22	26.84	32.76	37.42	44.27	49.22	54.07	58.20
D	6.00	8.00	10.00	12.00	15.00	18.00	23.00	30.00	36.00	45.00
C%	60.90	64.64	69.94	77.01	87.77	94.39	99.01	100.00	100.00	100.00
D	56.00	70.00	90.00	110.00	135.00	165.00	210.00	260.00	320.00	400.00
C%	100.00	100.00	100.00	100.00	100.00	100.00	100.00	100.00	100.00	100.00



C:\GRANULOG\OP36006.MES

CILAS 920 N.353

COMPAGNIE INDUSTRIELLE DES LASERS

H \* 4.23